

Softlin Komfort – more than a chair

The Softlin Komfort chair represents an optimal solution of sitting furniture for people with limited mobility. Its use is in hospitals as well as geriatric houses. Thanks to its wide modularity its possible to find such a variant of Softlin Chair, which will maximally meet your specific demands. In addition-ally your choice is extended by a variation of different accessories.


- HIGH FUNCTIONALITY
- MAXIMUM COMFORT
- ENJOYABLE DESIGN
- SAFETY
- EASY MANIPULATION
- EASY MAINTENANCE



Softlin Komfort – accessories



Crutches holder for practical placement of crutches while sitting



Holder for newspapers and magazines



Flexible hand support



Crutches holder with a holder for easy clamping



Revolving shelf with enhanced edges



Removable ergonomic eating plate with enhanced edges



Comfortable and practical removable covers of the armrest



Hook on the back part of the construction for fixation and storage of the eating desk



Hand approachable pocket in the back part of the backrest



Stowaway footboard with stowage compartment

TECHNICAL PARAMETERS

Width	70 cm
Height (depending on the back rest position)	98 - 113 cm
Depth (depending on the adjustment)	72 - 132 cm
Height of the sit section	47 cm
Diameter of castors	7,5 cm
5th central castor	10 cm
Lifting capacity	120 kg