

SOFTLIN PLUS

A CHAIR THAT WILL SIT YOU



SOFTLIN plus – a game for your comfort





SB 10101 – low



SB 20111 – medium



SB 40111 – low double chair



SB 30111 – high (without the head rest – option)

The SOFTLIN plus chairs

- **Complement the Softlin furniture range of bedside tables, beds with wooden lining, etc.**
- **With their technical parameters, structural design and hygienic properties they are mainly designed for hospitals, old people's homes, social care institutes, boarding houses for old people, home care, etc.**

- Robust wooden oak varnished structure
- Increased strength and stability
- The seat and back made of formed plywood
- Lifted seat and arm rests for easier sitting down and standing up
- Easy maintenance
- Three heights of the back
- Two positions of the arm rest
- Washable, removable textile overlay of the seat
- Self-extinguishing foam rubber
- A possibility of upholstery with artificial leather



Eating plate positioned on the arm rests (option).



Central nut for easy disassembly of the seat and removal of the overlay.




The design of the chairs allows you to stack them even in the case of arm rest models.



Head rest with adjustable height (recommended with the high back only)



Removable arm rests of oak lamellas, they can be adjusted to two heights, high rigidity of the arm rests provides for sufficient support during standing up. The arm rests can be assembled additionally to all types of SOFTLIN plus chairs.



The seat is coated with a separating watertight layer (we recommend to order spare overlays).

SB1	SB2	SB3	SB4 (SB5)

The manufacturer reserves the right to perform minor changes of the design that will not influence the function or the quality of the chair.

LINET spol . s r . o . , Želevčice 5, 274 01 Slaný, Czech Republic
Tel.: +420- 314/576 111, 524 139, fax: +420- 314/522 668
E-mail: info@linet.cz, http://www.linet.cz