

ProphyCair

Low Air Loss mattress



Patient

- Effective pressure distribution
- High level of patient safety
- Suitable for bariatric patients



Staff

- Adaptable to meet patient and clinical needs
- Infection control friendly



Management

- Cost effective pressure ulcer prevention
- Solution for heavier patients



Low Air Loss System

The highly effective low air loss system ProphyCair provides excellent pressure redistribution, patient immersion and envelopment. These aspects can significantly help in the prevention and treatment of pressure ulcers. The integrated Microclimate feature controls the moisture level of the skin aiding in the maintenance of skin integrity. The ProphyCair system is available in several sizes to suit all patients up to bariatric level of 445 kg.



9 different pressure settings according to the patient's needs

- 3** working modes
- Low Air Loss (LAL) mode – automatically maintains the selected pressure.
 - Pulse mode – every half minute reduces the pressure by 20% for 5 seconds.
 - Max Inflate – firm and stable support for patient care.

Safety of patient care

Fowler boost function

Once turned on, the mattress automatically increases the pressure if the backrest is raised to an angle of more than 35°. This reduces the risk of the patient bottoming out.



Patient safety

Audio and visual alarms in case of system faults. Self canceling with audio mute option.

Power fail alarm

If the mains power supply is lost when the system is in use.



Low pressure alarm

If the hose is disconnected or there is a major air leak.



Low Air Loss System

The ProphyCair Low Air Loss system encourages best possible levels of patient immersion and envelopment for effective skin pressure redistribution. In addition, it creates Microclimate air flow to aid the management of moisture around the patient's skin.



High quality materials

Durable and easy to clean and Dartex cover:

- waterproof
- vapour permeable
- 360° zipper

Mattress versions

The ProphyCair mattress is made in 4 different dimensions, all suitable for bariatric patients.

VERSION	DIMENSIONS	WEIGHT CAPACITY
ProphyCair™ Lo	200 × 87 × 18 cm	445 kg
ProphyCair™ Lo+	200 × 99 × 18 cm	445 kg
ProphyCair™	200 × 87 × 20.5 cm	445 kg
ProphyCair™+	200 × 99 × 25 cm	445 kg

The weight of the mattress is approx. 7 kg (depends on the type).



Members of LINET Group

LINET spol. s r. o.

Želevčice 5 | 274 01 Slaný | Czech Republic

tel.: +420 312 576 400 | fax: +420 312 522 668 | e-mail: info@linet.com | www.linet.com

