

■ ■ wissner-  
■ ■ bosserhoff


Functional elegance in nursing care

wissner-bosserhoff GmbH  
Hauptstraße 4-6  
D-58739 Wickede (Ruhr)  
Phone: +49(0)2377/784-0  
Fax: +49(0)2377/784-163  
E-Mail: info@wi-bo.de  
www.wi-bo.de

wissner-bosserhoff | Art.-No. 97-0744 | Edition: 05/2012, Rev.0 | Printer: msp druck und medien | Slight deviations in color are possible. Technical changes reserved.

■ ■ wissner-  
■ ■ bosserhoff


Functional elegance in nursing care



carisma 300  
Bariatric nursing home bed

# carisma 300

The bariatric bed is designed to provide a comfortable and safe care situation for residents with a weight of up to 270 kilograms. The extra-wide patient surface area of 100 cm ensures that residents do not feel restricted. Nursing staff have sufficient space for relocation, mobilization and positioning of residents and are provided with additional support with reinforced drives for bed height, backrest and leg rest adjustments.




Safe working load  
up to 300 kg

Like all nursing home beds from wissner-bosserhoff, the carisma 300 bariatric offers many user friendly features and a comfortable design in addition to the large patient surface area and the safe working load. For detailed information, please refer to the „carisma“ product catalog.

## Complementary products and accessories



The belsedia „big“  
bariatric chair  
Mod.-No. 60-0195



HeavyMatt 160 mattress  
Height 14 cm  
Mod.-No. 50-0315



HeavyMatt 200 mattress  
Height 18 cm  
Mod.-No. 50-0316

## Your advantages at a glance:

- Investment security:
  - Maximum resident weight up to 270 kg
  - Integrated bed extension
  - Minimum accessory requirements
- Comfortable patient surface area width of 100 cm
- Optimal egress and working heights individually adjustable between 36-80 cm
- Individual mobilization with divided side guards
- Load relief of the nursing staff