

Bed Magna 500



Bariatric Bed

Magna 500 - comfort for bariatric patients

The Olympia Bariatric bed combines an exceptional functionality and a high carrying capacity to support the treatment of bariatric patients.

- IDEAL FOR VERY HEAVY PATIENTS
- EASY TRANSPORT
- LIFTING CAPACITY OF 500 KG

CARRIAGE

- Very stable high/low adjustment by means of a support arm system (lifting capacity of 500 kg).
- Trendelenburg and Reverse Trendelenburg adjustment up to ±17°.
- High quality heavy-duty Twin castors Ø 150 mm with a directional castor.
- The brake pedal is accessible in all bed positions.
- The open carriage improves the accessibility under the bed and the ease of cleaning.

OPTIONS

- Atmosphere sideguards.
- Nurse box (ACO) + extendible drawer.
- Illuminated hand switch.
- Nightlight under the bed.

MATTRESS PLATFORM

- Mattress platform with 4 sections, adjustable in width (91 - 111 - 131).
- Backrest and knee rest are electrically adjustable.
- The electric manipulation of the mattress platform section is equipped with a clamp safety.
- Manual quick release of the backrest.
- Battery with alarm signal.
- Mattress platform with solid laminate sections of 6 mm, equipped with cutouts for patients restraint straps.
- Backrest and footrest are equipped with
 2 integrated mattress holders.
- Telescopic bed extension up to +20 cm.
- Equipotential connection.
- Bed corners with integrated bumpers.
- Cases to install an IV rod or a lifting pole.Bed panel with a solid laminate filling
- panel, choice out of 3 different width (104-124-144).
- The side beams on both sides are equipped with 3 synthetic accessory hooks.

ELECTRIC BED MANIPULATIONS

- The electric functions can be activated by means of a hand switch or a nurse box (*).
- The hand switch is equipped with a preprogrammed sleeping function.
- A nurse box (*) makes it possible to switch off several functions selectively.

TECHNICAL PARAMETERS

Outer dimensions	214 cm
(adjustable width)	104 -124 -144 cm
Mattress platform	206 cm
(adjustable width)	91-111-13# cm
Mattress platform height adjustment	42-82 cm
Backrest tilt	
(with autoregression)	±70°
Knee rest tilt	±34°
Weight (with options)	209 kg
Safe working load	500 kg

